

ReSound Smart Fit cable guide

Programming interface overview


ITEM NUMBER	DESCRIPTION	1	2	3	4	5	6
1	20310800	Noahlink Wireless	CS53 flex strip & CS44 socket cable (R&L)	Programming adapter & CS44 socket cable (R&L)	ENZO Programming adapter & CS44 socket cable (R&L)		
2	9022 907 69019	CS44 – Red, Long					
3	9022 907 69029	CS44 – Blue, Long					
4	19533900	CS53 Flex Strip					
5	19902500	PRG ADAPTER 2					
6	19215000	Programming Kit, ENZO					
ReSound OMNIA™		RIE, BTE, Rechargeable Customs	CIC	-	-		
ReSound ONE™		RIE, BTE	-	-	-		
Customs by ReSound		Rechargeable Customs	CIC	-	-		
ReSound LiNX Quattro™		Customs, RIE, BTE	Customs, RIE 61-DRW, BTE 67-DW	RIE 62-DRW, BTE 77-DW, BTE 88-DW	-		
ReSound ENZO Q™		BTE	-	BTE 88-DWH	BTE 98-DW		
ReSound Key™		Customs, RIE, BTE	Customs, RIE 61-DRW, BTE 67-DW	RIE 62-DRW, BTE 77-DW, BTE 88-DWH	BTE 98-DW		

NOTE: Strongly recommend to use fresh batteries for non-rechargeable models.

Cable connection for non-wireless programming

Custom models

With the battery door ajar, insert the flex strip into the slot beside the hinge pin. The black dot on the flex strip should be facing away from the battery door. The flex strip is fully inserted when the black dot abuts the hearing instrument case. Close the battery door once the flex strip is inserted.



RIE 62-DRW, BTE 77-DW & BTE 88-DW/DWH

To connect the programming adaptor, hold it at a 45° angle and insert the small pins at the tip of the adaptor into the square cut-out on the underside of the hearing instrument. Once in place, move the adaptor in the direction of the battery compartment and then gently click it into place.



RIE 61-DRW

Locate the programming slot below the push-button. Open the battery door. With the black dot on the flex strip facing the push button on the hearing instrument, insert the flex strip into the slot and close the battery door.



BTE 98-DW

[1] Locate and push out the battery door pin using the pin removal tool. [2] Pull out the battery door and [3] add the ENZO programming adaptor (19215000) and secure it with the pin. [4] Insert a new size #675 battery and close the door. [5] Connect CS44 socket cable.



Mini BTE 67-DW

Locate the programming slot below the microphone port. With the black dot on the flex strip facing the push button on the hearing instrument, insert the flex strip into the slot, until the black dot abuts the hearing instrument case.



ReSound battery guide

ReSound GN



BATTERY	BEHIND-THE-EAR HEARING AIDS						
	RECEIVER-IN-EAR			TRADITIONAL BEHIND-THE-EAR			
	Mini RIE		RIE	Mini BTE	BTE & Power BTE		Super Power BTE
	Integrated Lithium-ion	#312	#13	#312	Integrated Lithium-ion	#13	#675
ReSound OMNIA™	RUx60-DRWC RUx61-DRWC	RUx61-DRW	RUx61-DRW	-	RUx77-DWC RUx88-DWC	-	-
ReSound ONE	RTx61-DRWC	RTx61-DRW	RTx62-DRW	-	RTx77-DWC RTx88-DWC RTx88-DWHC	-	-
ReSound LiNX Quattro	REx61-DRWC	REx61-DRW	REx62-DRW	REx67-DW	-	REx77-DW REx88-DW	-
ReSound ENZO Q	-	-	-	-	-	EQx88-DWH	EQx98-DW
ReSound Key	-	KEx61-DRW	KEx62-DRW	KEx67-DW	-	KEx77-DW KEx88-DWH	KEx98-DW



BATTERY	CUSTOM HEARING AIDS						
	RECHARGEABLE CUSTOM		MICROPHONE-IN-HELIX	TRADITIONAL CUSTOM			
	Half Shell/ITC*	Full Shell	MIH	CIC	ITC	ITE/Half Shell	ITE/Full Shell
	Integrated Lithium-ion	Integrated Lithium-ion	#13/312	#10A	#312	#13/312	#13/312
ReSound OMNIA	RUxITC-DWC	RUxITE-DWC	-	RUxCIC	-	-	-
Customs by ReSound	RTxITC-DWC	RTxITE-DWC	-	RTxCIC	-	-	-
ReSound LiNX Quattro	-	-	RExMIH-W	RExCIC-W	RExITC-DW	RExITE-DW	RExITE-DW
ReSound Key	-	-	-	KExCIC-W	KExITC-DW	KExITE-DW	KExITE-DW

*Devices will be built as small as possible based on the patient's ear anatomy.

To ensure successful programming, always use a fresh battery. Remove the protective foil from the battery and wait for one to two minutes before inserting the battery into the hearing aid to allow the battery to fully activate.

ReSound Canada

2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
pro.resound.com



MK605113 REV K 2023.01