ReSound GN

ReSound Assist

# Snabbguide

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
The screenshot shows the ReSound Smart Fit 1.6 software interface. The top bar includes the title 'ReSound Smart Fit 1.6' and a search field for 'Hearing Care Professional'. The main interface is divided into several sections:

- Patient Information:** Nungo Hoegh, First Time User.
- Instruments:** Right and Left ReSound LINX Quattro 9 with RE961-DRWC MP Receivers. Both are connected and muted.
- Media Player:** Collapse Sidebar button.
- Patient Experience Level:** First Time User.
- Initial Programs:** Program 1: All-Around, Program 2: Restaurant, Program 3: None, Program 4: None.
- Right Instrument:** ReSound LINX Quattro 9, RE961-DRWC, MP Receiver, Serial Number: 9879879878, Cable: Wireless, Battery: Rechargeable.
- Left Instrument:** ReSound LINX Quattro 9, RE961-DRWC, MP Receiver, Serial Number: 9879879879, Cable: Wireless, Battery: Rechargeable.
- Audiogram:** A graph showing hearing levels (HL) across frequencies (125, 250, 500, 1K, 2K, 4K, 8K Hz). The graph is titled 'Audiogram was created 31.1.20'.

A red circle highlights the 'Hearing Care Professional' search field in the top right corner.

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1. Klicka på ”Start Consent”

**Patient Consent Needed**  
The sharing of health information requires patient consent. Review the document with the patient and Start Consent.

**Cloud-enabled dispensers**  
**Remote fine tuning of hearing instrument**

Contact details of your hearing care professional:  
Hearing Care Professional  
GN RESOUND US - MISC SALES  
8001 BLOOMINGTON FWY  
BLOOMINGTON  
MN 55420  
US

Details of data protection officer (if any) at hearing care professional:  
(Hereinafter "Hearing Care Professional" or "HCP")

Your Hearing Care Professional has sold you a hearing instrument produced by GN Hearing A/S, which is a Danish company with company reg. no. 55082715 and registered address at Lautrupbjerg 7, 2750 Ballerup, Denmark (hereinafter "the Manufacturer")

To access all of the features included with your hearing instrument, you will need to activate GN Online Services. These services enable you to receive firmware updates, and enable your Hearing Care Professional to fine-tune your hearing device remotely at your request.

By activating the GN Online Services, you will be allowing your Hearing Care Professional and the Manufacturer to access and process information related to your hearing loss via these online services and the associated mobile app. This notice explains what information will be collected, and asks for your consent to process this information for the purposes stated below.

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**Patient Consent Needed**  
The sharing of health information requires patient consent. Please click Agree to continue.

**PART A CONSENT**

Where you provide Part A Consent, you consent to your Hearing Care Professional and/or the Manufacturer (as specified) receiving and processing the following personal information about you in the manner, and for the purposes, detailed below which **IS NECESSARY** in order for the GN Online Services to be provided to you:

The following data categories are processed by the Manufacturer and HCP:

- **Your name and email address** : for the purpose of enabling fine-tuning and logging your consent.
- **Hearing threshold, hearing threshold shift, and other details regarding your hearing loss (health information)**: for the purpose of enabling fine-tuning.
- **The length of time you use your hearing instrument** : for the purpose of enabling fine-tuning and optimizing your settings.
- **The hearing instrument's serial number, hardware identification number and software version** for the purposes of enabling fine-tuning and enabling firmware updates.
- **Fitting specifications and fitting type** : for the purpose of enabling fine-tuning and optimizing your settings.
- **Hearing instrument settings** : for the purpose of enabling fine-tuning and optimizing your settings.
- **Information about the sound environment when sending a request via the app** : for the purpose of enabling fine-tuning and optimizing your settings.

The following data categories are processed by the Manufacturer:

- **Device type information** : Generic device-specific information such as the type of mobile device and mobile operating system: for the purpose of optimizing your settings and providing customer services (questions, complaints, repair, etc).
- **Technical log information** , Technical information in server logs, e.g., details of how the apps are used, Internet protocol, device event information such as crashes, system activity, hardware settings, browser type, browser language, the date and time of your request

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**Patient Consent Summary**  
Please print and sign the document below. This should be in the Hearing Care Professional's records.

First Name: Nungo Last Name: Hoegh Patient Email: patient@email.com

**Cloud-enabled dispensers**  
**Remote fine tuning of hearing instrument**

Contact details of your hearing care professional:  
Hearing Care Professional  
GN RESOUND US - MISC SALES  
8001 BLOOMINGTON FWY  
BLOOMINGTON  
MN 55420  
US

Details of data protection officer (if any) at hearing care professional:  
(Hereinafter "Hearing Care Professional" or "HCP")

Your Hearing Care Professional has sold you a hearing instrument produced by GN Hearing A/S, which is a Danish company with company reg. no. 55082715 and registered address at Lautrupbjerg 7, 2750 Ballerup, Denmark (hereinafter "the Manufacturer")

To access all of the features included with your hearing instrument, you will need to activate GN Online Services. These services enable you to receive firmware updates, and enable your Hearing Care Professional to fine-tune your hearing device remotely at your request.

By activating the GN Online Services, you will be allowing your Hearing Care Professional and the Manufacturer to access and process information related to your hearing loss via these online services and the associated mobile app. This notice explains what information will be collected, and asks for your consent to process this information for the purposes stated below.

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


The image shows three sequential screenshots of the 'Request assistance' form:

- Screenshot 1:** 'WHEN IS IT OCCURRING?' section. The user has selected 'With specific sounds' and 'Loud sounds: clanging dishes or cutlery, doors slamming, flushing water'.
- Screenshot 2:** 'HOW SEVERE IS THE PROBLEM?' section. The user has selected 'Slightly annoying'. Below, 'WHICH SIDE IS THE PROBLEM ON?' shows 'Both' selected. 'WHICH PROGRAM(S)?' shows 'All programs' selected.
- Screenshot 3:** 'HOW SEVERE IS THE PROBLEM?' section. The user has selected 'Slightly annoying'. Below, 'WHICH SIDE IS THE PROBLEM ON?' shows 'Both' selected. 'WHICH PROGRAM(S)?' shows 'All-Around' and 'Restaurant' selected.

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The image shows two final screenshots of the 'Request assistance' form:

- Screenshot 1:** 'SUMMARY' screen. It displays the selected options: 'Too sharp, tinny, or shrill', 'With specific sounds: Loud sounds: clanging dishes or cutlery, doors slamming, flushing water', 'Slightly annoying', 'Ear side: Both', and 'Restaurant'. Below, the 'REQUEST TITLE' field contains 'Too sharp, tinny, or shrill' and the 'PERSONAL MESSAGE' field contains 'Too sharp in restaurants'. A red box highlights the 'Send request' button.
- Screenshot 2:** 'THANK YOU!' screen. It displays a confirmation message: 'I have received your request for assistance. I will get back to you within 2 working days.' followed by 'Kind regards, Your hearing care professional' and the date '15 January 2017'. A red box highlights the 'Close' button.

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The screenshot displays the 'Fitting' tab of the software, which is divided into four sub-sections: Gain Adjustments, Advanced Features, Environmental Optimizer, and Device Controls. The Gain Adjustments section contains two frequency response graphs. The left graph shows a shaded area representing the hearing aid's performance and several colored lines with diamond markers representing individual patient settings. The right graph shows a similar view with a different set of patient settings. Below the graphs is a 'Gain Level (%)' dropdown menu set to 100%. The Advanced Features section below contains two identical grids of frequency response data. Each grid has columns for frequencies (All, 250, 500, 750, 1K, 1.5K, 2K, 3K, 4K, 6K) and rows for different hearing aid settings (50, 65, 80, CR, MPO). Below these grids are 'Handles' (6, 9, 17) and 'Compression' (WDRC) settings. At the bottom of the fitting screen, there are buttons for 'Calibrate DFS', 'Pair Accessories', 'On Off View Patient Settings', and 'Save'. The 'On Off View Patient Settings' button is highlighted with a red box and labeled '4'. The 'Noise Tracker II' section at the very bottom shows a slider with settings 'Off', 'Mild', 'Moder...', 'Consid...', and 'Strong Per Env...'. The 'Strong Per Env...' setting is highlighted with a red box and labeled '2'.

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

20 Feb 2018

**In-Office Session - 11:57 AM** Rating: Not satisfied yet

Action	Action by	Date	Time
Rate my sound Not satisfied yet	J. Snow	2/20/2018	11:57 AM <a href="#">Actions</a>
In-Office Session Saved in cloud	C. Chan	2/13/2018	11:35 AM <a href="#">Actions</a>

20 Feb 2018

**In-Office Session - 11:57 AM** Rating: Not satisfied yet

Action	Action by	Date	Time
Rate my sound Not satisfied yet	J. Snow	2/20/2018	11:57 AM <a href="#">Actions</a>

**Rating information from patient**  
How satisfied are you with your sound settings?

Dissatisfied      Not satisfied yet      Satisfied

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GN ReSound A/S
GN Hearing Sverige AB
GN Hearing Finland Oy/Ab