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1. Picou EM. MarkeTrak 10 (MT10) Survey Results Demonstrate High Satisfaction with and Benefits from Hearing Aids. *Seminars in Hearing*. 2020; 41(1):21-36.
2. Groth J. Hearing aid directionality with binaural processing. *AudiologyOnline*. 2016 May. Available from [www.audiologyonline.com](http://www.audiologyonline.com).
3. Appleton J, König G. Improvement in speech intelligibility and subjective benefit with binaural beamformer technology. *Hearing Review*. 2014;21(10):40-2.
4. Picou EM, Aspell E, Ricketts TA. Potential benefits and limitations of three types of directional processing in hearing aids. *Ear and Hearing*. 2014 May 1;35(3):339-52.
5. Jespersen CT, Kirkwood B, Groth J. Effect of directional strategy on audibility of sounds in the environment for varying hearing loss severity. *Canadian Audiologist*. 2017;4(6). Available from: <http://canadianaudiologist.ca/issue/volume-4-issue-6-2017/directional-strategy-feature/>.
6. Best V, Roverud E, Mason CR, Kidd Jr G. Examination of a hybrid beamformer that preserves auditory spatial cues. *The Journal of the Acoustical Society of America*. 2017 Oct 12;142(4):EL369-74.
7. Brimijoin WO, Whitmer WM, McShefferty D, Akeroyd MA. The effect of hearing aid microphone mode on performance in an auditory orienting task. *Ear Hear*. 2014; 35(5): e204-e212.
8. Archer-Boyd AW, Holman JA, Brimijoin WO. The minimum monitoring signal-to-noise ratio for off-axis signals and its implications for directional hearing aids. *Hearing Research*. 2018 Jan 1;357:64-72.
9. Van den Bogaert T, Klasen TJ, Moonen M, Van Deun L, Wouters J. Horizontal localization with bilateral hearing aids: Without is better than with. *The Journal of the Acoustical Society of America*. 2006 Jan;119(1):515-26.
10. Groth J, Laureyns M, Piskosz M. Double-blind study indicates sound quality preference for surround sound processor. *Hearing Review*. 2010;17(3):36-41.
11. Goyette A, Crukley J, Galster J. The Effects of Varying Directional Bandwidth in Hearing Aid Users' Preference and Speech-in-Noise Performance. *American Journal of Audiology*. 2018 Mar 8;27(1):95-103.
12. Groth J, Laureyns M. Preserving localization in hearing instrument fittings. *The Hearing Journal*. 2011 Feb 1;64(2):34-8.
13. Møller K, Jespersen C. The Effect of Bandsplit Directionality on Speech Recognition and Noise Perception. *Hearing Review Products*. 2013 Jun:8-10.
14. Pavlovic CV. Band importance functions for audiological applications. *Ear and Hearing*. 1994 Feb;15(1):100-4.
15. Ricketts TA. Directional hearing aids. *Trends in Amplification*. 2001 Dec;5(4):139-76.
16. Bo Nielsen J, Dau T, Neher T. A Danish open-set speech corpus for competing-speech studies. *The Journal of the Acoustical Society of America*. 2014 Jan;135(1):407-20.
17. Walden B, Surr R, Cord M, Dylund O. Predicting hearing aid microphone preference in everyday listening. *J Am Acad Audiol*. 2004;15;365-96.
18. Walden B, Surr R, Cord M, Grant K, Summers V, Dittberner A. The robustness of hearing aid microphone preferences in everyday environments. *J Am Acad Audiol*. 2007;18;358-79.
19. Preves DA, Sammeth CA, Wynne MK. Field trial evaluations of a switched directional/omnidirectional In-The-Ear hearing instrument. *Journal of the American Academy of Audiology*. 1999 May 1;10(5):273-84.
20. Humes LE, Rogers SE, Main AK, Kinney DL. The acoustic environments in which older adults wear their hearing aids: insights from datalogging sound environment classification. *American Journal of Audiology*. 2018 Dec 6;27(4):594-603.
21. Rallapalli V, Anderson M, Kates J, Balmert L, Sirow L, Arehart K, Souza P. Quantifying the Range of Signal Modification in Clinically Fit Hearing Aids. *Ear and Hearing*. 2020 Mar 1;41(2):433-41.
22. Orton JF, Preves D. Localization as a function of hearing aid microphone placement. *Hearing Instruments*. 1979; 30(1); 18-21.
23. Westerman S, Topholm J. Comparing BTEs and ITEs for localizing speech. *Hearing Instruments*. 1985; 36(2); 20-24.
24. Udesen J, Piechowiak T, Gran F, Dittberner A. Degradation of spatial sound by the hearing aid. *Proceedings of ISAAR 2013: Auditory Plasticity – Listening with the Brain. 4th Symposium on Auditory and Audiological Research*. August 2013, Nyborg, Denmark. Dau T, Santurette S, Dalgaard JC, Tanejbaerg L, Andersen T, Poulsen T eds.
25. Kollmeier B, Peissig J, Hohmann V. Real-time multiband dynamic range compression and noise reduction for binaural hearing aids. *Journal of Rehabilitation Research and Development*. 1993; 30: 82-94.
26. Carette E, Van den Bogaert T, Laureyns M, Wouters J. Left-right and front-back spatial hearing with multiple directional microphone configurations in modern hearing aids. *J Am Acad Audiol* 2014;25(9):791-803.

27. Groth J. The technical proof for clearer, fuller and richer sound with ReSound LiNX Quattro. ReSound white paper. 2018.
28. Groth J. An innovative RIE receiver with microphone in the ear lets users “hear with their own ears”. ReSound white paper. 2020.
29. Cord MT, Surr RK, Walden BE, Dittberner AB. Ear asymmetries and asymmetric directional microphone hearing aid fittings. American Journal of Audiology. 2011.
30. Zurek PM. Binaural advantages and directional effects in speech intelligibility. In G. Studebaker & I. Hochberg (Eds.), Acoustical Factors Affecting Hearing Aid Performance. Boston: College-Hill, 1993.
31. Cord MT, Walden BE, Surr RK, Dittberner AB. Field evaluation of an asymmetric directional microphone fitting. J Am Acad Audiol. 2007;18:245-56.
32. Bentler RA, Egge JLM, Tubbs JL, Dittberner AB, Flamme GA. Quantification of directional benefit across different polar response patterns. J Am Acad Audiol. 2004;15:649-59.
33. Hornsby B. Effects of noise configuration and noise type on binaural benefit with asymmetric directional fittings. Seminar presented at: 155th Meeting of the Acoustical Society of America; June 30-July 4, 2008; Paris, France.
34. Kirkwood B, Jespersen CT. How asymmetric directional hearing aid fittings affect speech recognition. Canadian Audiologist;4(1). Available from: <https://www.canadianaudiologist.ca/issue/volume-4-issue-1-2017/asymmetric-speech-recognition-feature/>.
35. Drullman R, Festen JM, Plomp R. Effect of temporal envelope smearing on speech reception. The Journal of the Acoustical Society of America. 1994 Feb;95(2):1053-64.
36. Kluender KR, Coady JA, Kieft M. Sensitivity to change in perception of speech. Speech Communication. 2003 Aug 1;41(1):59-69.
37. Reinhart P, Zahorik P, Souza P. The interaction between reverberation and digital noise reduction in hearing aids: Acoustic and behavioral effects. The Journal of the Acoustical Society of America. 2017 May;141(5):3971.