

ReSound Smart Fit 1.15

Improvements for fitting efficiency

Provide your best hearing care, even remotely, fitting many needs and preferences efficiently with our optimized navigation flow and fitting process in ReSound Smart Fit™ 1.15.

What's new in ReSound Smart Fit 1.15?

This update supports your daily fittings in becoming more accurate and intuitive by improving:




Fitting Flow & Navigation

Verification

Personalization



Fitting Flow & Navigation

Instant hearing instrument detection and easy access to key actions enable you to provide a great first fit:

What's new?	What it is?	Benefit	Applicable to	How does it look?
Assign Instruments screen upgrades	<p>Re-organizing the "Assign Instrument" screen by further separating the hearing instruments assigned from the others that have been detected by the software</p> <p>Changing the handling of the "Continue" button by disabling it until you are properly set up for success in the next steps</p>	<p>Easy selection in the beginning of the fitting so you can clearly see what you assigned.</p> <p>Restrictions around "Continue" button enable you to make the right choice before continuing the fitting.</p>	All ReSound hearing instruments in ReSound Smart Fit	
Push-to-Assign	<p>Identify and assign hearing instruments to left and right side with a single push on the program button</p>	<p>Assigning hearing instruments in a single step saves you the number of steps to start connecting to the hearing instruments</p> <p>An alternative to using beeps and clicking to identify and assign the hearing instruments</p>	ReSound OMNIA	
Updated navigation	<p>Easier accessibility to fine tuning of sound processing features and device controls.</p> <p>Promoting previously hard-to-find features e.g. In-Situ Fitting and TSG.</p> <p>Relabeling screens for a more intuitive access and removing redundancies.</p>	<p>Freeing up time for hearing care and reducing time spent on navigating the software</p>	All ReSound hearing instruments in ReSound Smart Fit	



Verification

Give the best sound experience from the start with more helpful indicators to support your daily fittings:

What's new?	What it is?	Benefit	Applicable to	How does it look?
AutoREM	Physical Properties overview has been added to the first step of the AutoREM flow Probe Tube placement indicator available to reassure all is in place when your patient is all set	Supports real-ear measurement use	All ReSound hearing instruments in ReSound Smart Fit	
Fit to Target	If you use the NAL NL2 fitting rule, Smart Fit 1.1.5 sets Speech Weighted Noise as initial gains	Easier to match targets and greater fitting accuracy for when using NAL-NL2 as fitting rule	ReSound OMNIA	

Personalization

Intuitive, easy-to-use enhancements to support personalized care for your patients:

What's new?	What it is?	Benefit	Applicable to	How does it look?
Automated Preferences	Smart Fit 1.1.5 keeps count on the number and share of your sessions and changes automatically certain settings to your preferences. When a certain number of sessions with redundant steps have been reached, Smart Fit will inform you about the update of your fitting preferences.	Automated preference update to remove redundant steps performed in the sessions	All ReSound hearing instruments in ReSound Smart Fit	
Receiver Recommendation	Receiver recommendation based on your patient's audiogram available in Smart Launcher screen	Helps you to select the appropriate receiver for your patient	Helps you to select the appropriate receiver for your patient	

If you have questions about ReSound Smart Fit and the new updates, please contact Technical Support.

Manufacturer according to FDA:

GN ReSound North America
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-248-4327
pro.resound.com

ReSound Government Services
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
gs.resound.com

Manufacturer according to Health Canada:

ReSound Canada
2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
pro.resound.com

