

People with hearing loss are at the heart of everything we do because we know they're placing their trust in us. This is why we've been pioneering innovative new technologies, with consistent quality, for over 150 years to help people feel more involved, connected and in control, transforming their lives through the power of sound. Available exclusively through audiologists and hearing care professionals.

ReSound is part of the GN Group – pioneering great sound from world-leading ReSound hearing aids to Jabra office headsets and sports headphones. Founded in 1869, employing over 5,000 people, and listed on NASDAQ OMX Copenhagen, GN makes life sound better.



Manufacturer according to FDA:

GN ReSound North America
8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-882-3636
future-resound.com

Manufacturer according to Health Canada:

ReSound Canada
2 East Beaver Creek Road, Building 3
Richmond Hill, ON L4B 2N3
Canada
1-888-737-6863
future-resound.com

© 2020 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Google Play and the Google Play logo are trademarks of Google Inc.



MK605223 Rev B 2020.04



ReSound Preza™

Caring for your rechargeable hearing aids





GN Making Life Sound Better

Charging

ReSound Preza features a completely sealed lithium-ion battery. With its small, discreet size and long-lasting battery, you get power that lasts more than a full day – without having to change a single battery. In addition, the charging case’s built in battery holds three full charges of power, giving you the freedom to charge anywhere.

Upon insertion or removal of hearing instruments, the current battery status is indicated for 5-10 seconds



Charging case battery status

HOW TO CHARGE RESOUND PREZA

1. Make sure the hearing instruments are dry and clean before putting them in the charging case.
2. When inserted in the charging case correctly, the five LED lights on the front of the charging case will indicate the current battery level. The LED on the hearing instruments slowly pulsates throughout the charging process. When fully charged, the LED will be solid green until removed from the charger.
3. Gently pull the hearing instrument by the casing to remove from the charging case. The hearing instruments are automatically powered on when removed.

CHARGING TIME	BATTERY LIFE
3 hours =	30 hours 24 hours with 50% streaming time
1 hour =	16 hours
30 minutes =	8 hours

HOW TO CHARGE THE PORTABLE CHARGING CASE

When charging your portable charging case for the first time, charge it for at least three hours, even if the battery status LED lights indicate that the battery is fully charged.

1. Plug the supplied power adaptor into a power outlet and connect the charging cable to the back of the charging case.

A close-up photograph of a dark-colored hearing aid device. The text "ReSound GN" is printed in a light, sans-serif font on the curved surface of the device. The background is white, and the device is the primary focus of the top and bottom portions of the page.

ReSound GN

CHARGING CONSIDERATIONS

- It is recommended to charge your hearing instruments every night and it is safe to leave the charging case connected during that time. The batteries in both the hearing instruments and portable charging case cannot be overcharged.
- Charging should take place between 10 °C (50 °F) and 40 °C (104 °F). Your hearing aids may not charge outside of this range, however, they can be worn outside of these limits. Keep your charger out of direct sunlight or very warm temperatures while charging your hearing aids, just as you would with your mobile phone or other rechargeable devices.
- If the charging case runs out of power during the charging process, the hearing instruments will power off to preserve battery. When you remove the hearing instruments from the charging case, they will not power on automatically. To power on manually, press the push button for 5 seconds and the LED will light up.
- The hearing instruments will charge whether the charging case is open or closed.

Care and maintenance

Maintaining your hearing instruments through daily cleanings and regular service visits is important. Proper care helps retain optimum hearing conditions, extends the life of your hearing instruments and ensures good hygiene.

Tips:

- Have the Costco Hearing Aid Center professionals show you how to properly clean your hearing instruments.
- Purchase hearing instrument cleaning tools through your Costco Hearing Aid Center.
- Hold the hearing instrument over a soft surface, in case it is dropped during cleaning.
- Never use water, alcohol, solvents or cleaning fluids on the hearing instruments.

DAILY

- Remove your hearing instruments when applying cosmetics, perfume, aftershave, hair spray and suntan lotion as they may get into the hearing instrument and cause damage.
- Clean your hearing instruments using a soft cloth, tissue or proper cleaning tools to remove earwax, dust or moisture.
- Store your hearing instruments in the charging case when not in use.

REGULARLY

- Use a damp cloth to clean the receiver tube and receiver dome.
- Remove dust or debris from the charging case inserts and receiver bay with clean cotton bud. Clean the outside surface with a soft cloth.
- Contact your Costco Hearing Aid professional for replacement components when they become stiff, brittle or discolored.

CARE AND MAINTENANCE CONSIDERATIONS

- If immersed in water, your hearing instruments should be removed immediately. Dry them with a clean towel and place them in a dry place. Do not place the hearing instruments in the charging case when wet.
- Hearing instruments should be professionally cleaned regularly as recommended by your Costco Hearing Aid professional to improve the life of your hearing instruments.
- If necessary to use heat dryers to dry your hearing aids, please ensure that your hearing aids are cooled off to room temperature prior to drying. The temperature of a dryer should not exceed 40 °C (104 °F) at any time during the drying cycle. Exposure to 40 °C (104 °F) temperatures and above can shorten the life and degrade the lithium-ion battery.

