



ReSound cable guide

ReSound Aventa® 3 		 		 			
 <p>Truly wireless fitting (no cable required)</p>		 <p>3-pin CS53 flex strip and 4-pin CS44 socket cable*</p>		  <p>Programming Adaptor I or Programming Adaptor II and 4-pin CS44 socket cable</p>		 <p>ENZO Programming Adaptor 19215000</p>	
Fresh battery		Fresh battery		Fresh battery			
ReSound LiNX2™	Wireless custom, RIE, BTE	Custom, RIE: 61, Mini BTE: 67		RIE: 62, BTE			
ReSound ENZO2™	Wireless BTE					BTE: 88	BTE: 98
ReSound Enya™	Wireless custom, RIE, BTE	Custom, RIE: 61		RIE: 62, BTE			
ReSound LiNX™	Wireless BTE, RIE: 61	RIE: 61		BTE			
ReSound ENZO™	Wireless BTE						BTE: 98
ReSound Verso™	Wireless custom, RIE, BTE	Wireless custom, RIE, BTE		RIE: 62, BTE			
ReSound Alera®							
ReSound Vea™		Custom		RIE: 62, BTE: 77, 88		BTE: 60, 70, 80	
ReSound Magna™						BTE	

1 Custom models

With the battery door ajar, insert the flex strip into the slot beside the hinge pin. The black dot on the flex strip should be facing away from the battery door. The flex strip is fully inserted when the black dot abuts the hearing instrument case. Close the battery door once the flex strip is inserted.

NOTE: ReSound LiNX2 and ReSound Enya customs require CS53 (19533900) flex strip for cabled programming.

2 RIE: LS61 and LN61

Locate the programming slot below the push-button. Open the battery door. With the black dot on the flex strip facing the push button on the hearing instrument, insert the flex strip into the slot and close the battery door.

3 Mini BTE: LS67





Locate the programming slot below the microphone port. With the black dot on the flex strip facing the push button on the hearing instrument, insert the flex strip into the slot, until the black dot abuts the hearing instrument case.

4 5 RIE: 62, and BTE models

To connect the programming adaptor, hold it at a 45° angle and insert the small pins at the tip of the adaptor into the square cut-out on the underside of the hearing instrument. Once in place, move the adaptor in the direction of the battery compartment and then gently click it into place.

6 BTE:EN/EO98

Remove the existing battery door and replace with the DA/battery door (19215000).

AVENTA 2.9x		 <p>4-pin CS44 socket cable</p>		 <p>3-pin CS53 flex strip and 4-pin CS44 socket cable*</p>		 <p>Programming Adaptor I and 4-pin CS44 socket cable</p>		 <p>Battery pill and 4-pin CS44 socket cable</p>	
Non-wireless programming		NO battery	Fresh battery	NO battery	NO battery	NO battery	NO battery	NO battery	NO battery
ReSound Live®			CUSTOM		BTE, RIE				
dot² by ReSound®			RIE						
dot by ReSound®									
ReSound Essence®	BTE		CUSTOM						
ReSound Sparx®	BTE								
be by ReSound®			CUSTOM					IOT™	
ReSound Azure®									
ReSound Xplore®	BTE: 60, 70, 80		CUSTOM		Silhouette BTE: 61, 71, 81				
ReSound Ziga®									

Reminder: To ensure successful programming, be sure to note whether a battery is required. When a battery is required, it is important to use a fresh battery while programming.

* The CS44 cable and CS53 flex strip combination can be used interchangeably with the CS63 cable and CS63 flex strip combination.

ReSound battery guide



Behind-the-ear instruments

Product	Receiver-in-ear				Standard				
	Micro RIE	Mini RIE	RIE		Mini BTE	Transitional BTE	BTE	Power BTE	Super Power BTE
	#10A	#312	#312	#13	#312	#13	#13	#13	#675
ReSound LiNX ² ™		LSx61		LSx62	LSx67		LSx77	LSx88	
ReSound ENZO ² ™								ENx88	ENx98
ReSound Enya™			EY362		EYx67		EYx77	EYx88	
ReSound LiNX™		LNx61					LNx77	LNx88	
ReSound ENZO™									EOx98
ReSound Verso™	VOx60 VOTx60	VOx61 VOTx61	VOx62 VOTx62		VOx67 VOTx67		VOx77 VOTx77	VOx88 VOTx88	
ReSound Alera®	ALx60 ALTx60	ALx61 ALTx61	ALx61 ALTx61		ALx67 ALTx67	ALx77/87 ALTx77/87			
ReSound Vea™			VE362		VEx60		VEx70 VE377	VEx70 VE377	
ReSound Magna™									MG490

Custom instruments

Product	Microphone-in-helix			Traditional custom					
	MIH			IIC	CIC	Mini canal	ITC	ITE half shell	ITE full shell
	#10A	#312	#13	#10A	#10A	#10A	#312	#312	#13
ReSound LiNX ² ™	LSxMIH-S	LSxMIH	LSxMIH	LSxIIC	LSxCIC		LSxITC	LSxITE	LSxITE
ReSound Enya™					EYxCIC		EYxITC	EYxITE	EYxITE
ReSound Verso™	VOx10-M VOTx30-M	VOx30-M VOTx30-M	VOx30-M VOTx30-M	VOx10-C VOT910-C	VOx10 VOT910		VOx30 VOT930	VOx50 VOT950	VOx50 VOT950
ReSound Alera®	ALx10-M	ALx30-M ALTx30-M	ALx50-M	ALx10-C	ALx10	ALx20	ALx30	ALx40	ALx50
ReSound Vea™					VEx10		VEx30	VEx50HS	VEx50

Manufacturer according to FDA:

ReSound
North America
8001 Bloomington Freeway
Bloomington, MN 55420
1-800-248-4327
resoundpro.com

ReSound
Government Services
8001 Bloomington Freeway
Bloomington, MN 55420
1-800-392-9932
resound.com/governmentservices

Manufacturer according to Health Canada:

ReSound
Canada
303 Supertest Road
Toronto, Ontario M3J 2M4
1-888-737-6863
resoundpro.com

