


# More ways to connect

Made for iPhone | iPad | iPod

Accessed through the ReSound Smart 3D™ app, ReSound Assist Live Assistance gives you a convenient way to connect with your audiologist through live, video appointments.

After your first fit at the hearing clinic, you'll enjoy continuous support and all the guidance you need in the most convenient and comfortable way for you and your audiologist. It's all designed to help you get the most out of your ReSound hearing aids.




## What is it?

ReSound Assist Live Assistance lets you supplement clinic visits with convenient, real-time video appointments. That means you can enjoy personalized, flexible care from the comfort of your own home. It's as easy as answering a video call on your phone.

## How does it work?

1. Book a video appointment with your audiologist. At the appointment time, your audiologist will initiate the call.
2. Answer the video call received through the ReSound Smart 3D app.
3. Your audiologist connects to your hearing aids and provides you with live adjustments, if needed.
4. You evaluate any new settings made with your audiologist in real-time during the appointment. When the call is complete, your hearing aids will automatically reboot and new settings will be saved.



Find out more about ReSound Assist Live Assistance at [resound.com/assist](https://resound.com/assist).

## Manufacturer according to FDA:

### ReSound Government Services

8001 E Bloomington Freeway  
 Bloomington, MN 55420  
 USA  
 1-800-392-9932  
[resound.com/veterans](https://resound.com/veterans)

