

ReSound OMNIA

So incredible we almost couldn't tell you – but now we can

Now your clients can hear their best in noise with an almost unbelievable 150% improvement in speech understanding.^{*} Powered by the industry's only Organic Hearing[™].

> *4dB SNR improvement to speech understanding compared to legacy products. Jespersen et al (2022)

Wow your clients from the start

Combining the strength of humans and machine with our cutting-edge technology and leading audiological expertise, we developed our most intelligent and automated steering and beamforming system yet. This gives users superior hearing clarity and understanding in all sound environments, so they can love the everyday details in life again.

OUR TRILOGY OF TECH

- How we achieved this audiological feat:
- 1. Widening our sound environment mapping (360 All-Around)
- 2. Speeding up our processing and steering for faster
- 3. Narrowing the beamforming to give even more precision for our best speech clarity in any environment

All working at the same time across any sound environment, enabling your clients to automatically tune into what they want to listen to and not have the hearing aids choose for them.




A SUCCESSFUL FIRST FIT, NATURALLY

Enjoy more time for your clients and work even more efficiently with our optimised fitting software and tools, which help you give your best hearing care, even remotely, and fit many needs and preferences.

Further improve your fitting success for all ReSound products by 13% with ReSound Smart Fit[™] 1.15, the very latest version of our fitting software.

Powered by Organic Hearing



SOUNDS NATURAL

Emulating the natural hearing process and the natural ways we listen as closely as possible, helping the brain to hear as nature intended.




Our best speech understanding in noise*1



Hearing clarity all around¹



Natural own voice perception²



FEELS NATURAL

Adapts organically to the person for the best wearing comfort all day long, with ease of use.



Optimal comfort on the ear



Check My Fit for easier onboarding

All-weatherproof design for less maintenance



CONNECTS YOU NATURALLY Inspired by the way we intuitively interact with people and devices in our daily lives.



-

Crystal clear phone calls

Hands-free calls for iPhone and iPad*

Direct streaming from your favourite devices

Compared to our legacy technology

- Jespersen et al (2022) Groth et al (2022)


* ReSound OMNIA hands-free calls are compatible with iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), and iPad mini (6th generation) or later, with software updates iOS 15.3 and iPadOS 15.3 or later.



IMPROVE ONBOARDING & USAGE

Check My Fit helps users to wear their hearing aids properly and get the most out of ReSound OMNIA™ each day, which can lead to fewer service visits. Available in the ReSound Smart 3D[™] app.





STYLISH COLOURS

Grey

Get a beautiful finish on an ergonomic design that raindrops just slide off.





Price point & key feature overview

	Feature Name	RU9	RU7	RU5
Comfort	Noise Tracker II	5 settings	3 settings	2 settings
	Wind Guard	3 settings	2 settings	1 setting
	Impulse Noise Reduction	3 settings	3 settings	1 setting
	Expansion	3 settings	2 settings	1 setting
	Environmental Optimizer	Volume and PNR	Volume only	Volume only
Speech Understanding	360 All-Around	•		
	Front Focus	•		
	All Access Directionality		•	
	Ultra Focus		•	
	M&RIE	•	•	•
	Adjustable Directional Mix	•		
	AutoScope Adaptive Directionality	•		
	MultiScope Adaptive Directionality		•	
	Spatial Sense		•	
	Binaural Directionality III & Spatial Sense			•
	Soft Switching	•	•	•
	Adaptive Directionality	•	•	•
	Fixed Directionality	•	•	•
Feedback Management	DFS Ultra III	•	•	•
	Music Mode	•	•	•
	Whistle Control	•	•	•
	Auto DFS	•	•	•
Convenience	Ear-to-ear synchronization	•	•	•
	Hands-free calls for iPhone & iPad	•	•	•
	Made for iPhone, iPad and iPod touch (MFi)	•	•	•
	Android streaming (ASHA)	•	•	•
	Wireless Accessories	•	•	•
Fitting Flexibility	Max Gain Handles	17	14	12 (no 8K handle)
	Max Number of Programs	4	4	4
	ReSound Assist	•	•	•
	ReSound Assist Live	•	•	•
	Low Frequency Boost	3 settings	2 settings	1 setting
	Amplification strategy (WDRC-Semi-Linear-Linear)	3 settings	3 settings	2 settings

READY TO HEAR THE DIFFERENCE?

Try ReSound OMNIA for yourself to experience the benefits and see why it's a natural fit for your portfolio and your clients. Let's create wow moments with the everyday details in life your clients love.

Visit pro.resound.com to find out more

Worldwide Headquarters GN ReSound A/S DK-2750 Ballerup Denmark Tel.: +45 4575 1111 resound.com

CVR no. 55082715

United Kingdom GN Hearing UK Ltd. Unit 13 Talisman Business Centre Bicester Oxon OX26 6HR Tel: +44 1869 352800 resound.com

Australia

GN Hearing Australia Pty Ltd Gate C, 19-25 Khartoum Road Macquarie Technology Park Macquarie Park NSW 2113 Tel.: (free) 1800 658 955 resound.com

New Zealand

GN Hearing New Zealand Limited Ground Floor, North Entrance 4 Fred Thomas Drive Takapuna, Auckland, 0622 Tel.: (free) 0800 900 126 resound.com

Singapore

GN Hearing Pte. Ltd. 456 Alexandra Road, #22-01, Singapore 119962 Tel: +65 6320 9388 resound.com

M102892GB-22.09-Rev.C

© 2022 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple and the Apple logo, iPhone, iPad and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC