




ReSound LiNX 3D™



The highest classification accuracy, compared to 5 leading hearing aid brands

ReSound LiNX 3D hearing aids are up to

**50%**


better at identifying speech in various environments\*

■ ReSound LiNX 3D    □ Other brands

So your clients can

hear 80% more surrounding sounds\*\*

A lower threshold for picking up speech that comes from beside and behind the user, resulting in a more natural listening experience



When speech is coming from the left side (at 90 degrees) **ReSound LiNX 3D** has the lowest threshold for picking up speech in a noisy background



recognize **40%** more speech in noise\*\*

Your clients can hear more of what's in front of them – without losing access to speech and sound coming from other directions



Sources: \*Groth (2016) \*\*Jespersen et al. (2016)

# Binaural Directionality III and Spatial Sense



Supports the brain's ability to separate and focus on specific sounds – without losing touch with the environment

ReSound LiNX 3D delivers all the benefits of directional technology while preserving the ability to hear surrounding sounds



Specialized algorithms help the brain detect where sounds are coming from, providing your clients with a more natural sense of sound