



ReSound ONE™, ReSound LiNX Quattro™
and ReSound ENZO Q™

Hear like no other

GN Making Life Sound Better

Made for
iPhone | iPad | iPod

Powered
android

A close-up, profile view of a woman with blonde hair, holding a seashell to her ear. She is looking towards the right. The background is a soft, out-of-focus light blue and white, suggesting an outdoor setting like a beach. The lighting is warm and natural, highlighting the texture of her skin and the shell.

Inspired by Organic Hearing

The ReSound Organic Hearing philosophy enables you to connect to the world around you in the most intuitive and natural way. We achieve this by developing solutions that work with your individual ear anatomy to more closely mimic how sounds in the environment are naturally collected and delivered to the brain.

This philosophy drives a history of innovations that have given hearing aid wearers the ability to instinctively select the sounds that are unique to your world.

Hear like no other with ReSound ONE

Just like a fingerprint, each person's hearing is ONE of a kind. This is because nature has designed the shape of your ear to perfectly capture the sounds around you.

ReSound ONE™ with M&RIE (Microphone & Receiver-in-Ear) delivers a new class* of hearing aid with an additional microphone inside your ear, allowing the unique shape of your ear to collect sounds, just as nature intended.

With convenient direct audio streaming, access to remote hearing care services, and industry-leading battery life, ReSound ONE truly keeps you connected to people and technology you need to become ONE with your world.



ReSound ONE provides you with a more complete and individualized sound experience with:

A MORE NATURAL SOUND AND IMPROVED SENSE OF SPACE


ReSound ONE with M&RIE is the first hearing aid* with three microphones, including one in the ear, to use your ear's natural ability to collect the information your brain needs to make sense of sound.

EASIER-TO-FOLLOW CONVERSATIONS IN DYNAMIC ENVIRONMENTS

All Access Directionality gives you access to all the sounds around you, putting you - and not your hearing aids - in charge of which sounds you want to focus on.

BEST 1:1 SPEECH UNDERSTANDING

Like an invisible beam of sound detection, Ultra Focus lets you decide when to concentrate on the speech in front of you while minimizing background noise - even in the toughest environments.



*Microphone in ear and two standard directional microphones, directionality features, wireless audio streaming.

Unmatched connectivity to adapt to your needs

ReSound ONE adapts to your lifestyle by connecting you to the world around you with direct audio streaming from both iOS® and Android™ devices. So your favorite audio book, music or even GPS directions can stream directly to your hearing aids without an intermediary device.

One-tap sound adjustments on your ReSound Smart 3D™ app and a versatile range of wireless accessories give you even more options to connect to the sounds you enjoy every day, like your favorite TV program.

Why do users love ReSound ONE with M&RIE?

- Provides the most preferred, natural sound quality¹
- Helps you better localize sound from the start and that ability to localize sound improves over time²
- Delivers better wind noise protection so being outdoors is more enjoyable³
- Provides effortless, non-streamed phone calls⁴

1. With mild hearing loss compared to pinna restoration and omnidirectional. Groth, 2020.

2. Compared to traditional microphone placement. Jespersen et al, 2021.

3. Compared to traditional microphone placement. Andersen et al, 2021.

4. Hartenstein & Quilter, 2021. Subjective experience of telephone listening with M&RIE.

Industry-leading rechargeability

ReSound ONE offers the best rechargeable battery life in the industry with up to 30 hours* of battery life or up to 25 hours with the ability to stream as much as you want.

*Expected battery life dependent on active features, the use of wireless accessories, hearing loss, device age and sound environment.



Designed for ultimate comfort, discretion and performance

Built for ultimate comfort and adaptability, the modern styling of ReSound ONE has just one touchpoint on the back of the ear and sits lower for greater discretion.

Available in RIE and now BTE styles, ReSound ONE comes in a variety of colors to blend with hair and skin tones, so you can wear ReSound ONE with confidence.

BTE STYLES



BTE 77
Rechargeable solution



Power BTE 88
Rechargeable solution



High Power BTE 88
Rechargeable solution

RIE STYLES



Mini Receiver-in-Ear
Rechargeable solution



Mini
Receiver-in-Ear



Standard
Receiver-in-Ear

RIE and BTE COLORS



Sparkling
Silver

Champagne

Gold

Bronze

Warm
Grey

Graphite

Espresso

Deep
Black



BTE Desktop Charger



Getting started with hands-free calling

ReSound has made it more convenient than ever for patients to stay in touch with the people in their lives. Now patients can enjoy hands-free phone and FaceTime calls for iPhone and iPad*, using state-of-the-art streaming technology.

Using ReSound ONE™ as a wireless headset, your patients can effortlessly take calls and talk without needing to hold the iPhone or iPad. They just tap the button on their hearing aids to answer a call and get the call streamed directly to their ears. With low battery consumption, they have the power to keep on talking for a full day.

*iPhone 11 or later, iPad Pro 12.9-inch (5th generation), iPad Pro 11-inch (3rd generation), iPad Air (4th generation), and iPad mini (6th generation), with software updates iOS 15.4 and iPadOS 15.4 or later.

ReSound LiNX Quattro™

A hearing experience beyond compare

Our latest complete family of premium hearing aids is beyond compare with richer sound quality, direct streaming to both Apple® and Android devices and intuitive rechargeability - all in one product. Combined with remote support and personalization with the ReSound Smart 3D app and ReSound Assist Live Assistance, you get an unrivalled hearing experience.



Rechargeable
RIE 61

Mini RIE 61

RIE 62



Mini BTE 67

BTE 77

Power BTE 88



CIC-W

MIH

ITC

ITE-HS

ITE-FS

ReSound LiNX
Quattro 61 Portable
Charging Case



Brilliant sound experience

CLEARER SOUNDS

ReSound LiNX Quattro utilizes the full range of the microphone allowing sound to be processed in a cleaner way - free of distortion.

+


FULLER SOUNDS

ReSound LiNX Quattro provides you with a larger range of sounds and access to nuances in tones that were previously unavailable.

=

RICHER EXPERIENCE

With improved clarity and access to more sounds, coupled with a stronger wireless signal, your favorite environments become richer with ReSound LiNX Quattro.



ReSound LiNX Quattro is preferred 95% of the time for music listening, compared to other premium hearing aids.*

- ReSound LiNX Quattro
- Other premium MFi hearing aids

*Jespersen et al, 2018



Intuitive rechargeability

ReSound LiNX Quattro offers intuitive rechargeability with a lithium-ion battery that is completely sealed for protection. With its small, discreet size and long-lasting battery, you get power that lasts more than a full day – without having to fuss with any batteries.

FULL CONFIDENCE WITH A COMPLETE CHARGE

At night, simply place your hearing aids in the easy-to-use charging case to prepare for the next day. If you do forget to charge, the case's built in battery holds three full charges of power, giving you the freedom to charge anywhere.

CHARGING TIME	HOURS OF BATTERY LIFE
3 hours =	30 with 50% streaming time 24
1 hour =	16
30 minutes =	8

ReSound ENZO Q™

Power plus performance

ReSound ENZO Q, our premium hearing aid family for severe-to-profound hearing loss, provides more than just the power you need. It brings together the connectivity, support and personalization you want for the best hearing experience in any environment.

With ReSound ENZO Q you get:

60%

More clarity of the sounds around you*

60%

Better speech understanding in noise**

*compared to binaural beamforming in hearing aids for severe-to-profound hearing loss (*Effect of directional strategy on audibility of sounds in the environment for varying hearing loss severity*. Jespersen C, Kirkwood B, Groth J., 2017).

**with Binaural Directionality III compared to omni directionality (*Binaural Directionality III compared to omni directional, data on file*)



Crystal clear,
comfortable,
high-quality
sound

Severe-to-profound hearing loss requires a complete hearing solution that is more than powerful. ReSound ENZO Q is created with the most advanced technology so you can enjoy clear and comfortable sound, extensive connectivity options, and a 20% longer-lasting battery.

EXPERIENCE THE POWER OF CLARITY

Clarity of sound is important in any environment. Personalize your hearing experience to any situation and know that ReSound ENZO Q technology is working to eliminate the annoying whistling or static noise even when you turn up the volume for comfort and clarity.

MORE OPTIONS FOR CONNECTIVITY AND DIRECT STREAMING

ReSound ENZO Q offers a complete ecosystem of connectivity which includes direct streaming capabilities from both iOS and Android devices.* Stream your favorite music, audio book, talk on the phone or even get GPS directions streamed to both of your hearing aids. It's as simple as using wireless stereo headphones.



High Power
BTE 88



Super Power
BTE 98

*Full streaming compatibility requires minimum Android v10 and Bluetooth 5.0 and that functionality has been enabled by the phone manufacturer and service carrier. ReSound cannot guarantee full compatibility with all Android devices.

ReSound Smart 3D™ app


Support and Personalization

Tailor your hearing experience to get the most out of your hearing aids with the ReSound Smart 3D app.

ReSound hearing aids automatically adjust to your environment, but you can personalize the details of your hearing experience even further.

With ReSound Smart 3D app, you will have access to:

- Sound Enhancements in a variety of environments
- Save Favorites and GeoTag Locations
- Find My Hearing Aids functionality
- Hearing aid advice and tips for sound optimization
- Rechargeable battery status updates and notifications



Sound controls available on Apple Watch®



More connections for better care

The complete ReSound portfolio of hearing aids includes access to ReSound Assist, so you can receive hearing care in the most convenient way.

FACE-TO-FACE ONLINE HEARING CARE

ReSound Assist Live Assistance makes it possible to have face-to-face fitting sessions from the comfort of your home. Simply answer the incoming video call from your hearing care professional and get real-time fine-tuning adjustments.

Always discreet and reliable

Every component of ReSound products, inside and out, is fully coated with iSolate nanotech to stay sealed and protected against dirt, moisture, dust or earwax. So you can step out into the rain or onto the field with total confidence.

ReSound RIE and BTE styles



ReSound Custom styles



Wireless accessories

ReSound has an entire portfolio of wireless accessories that work together as an intelligent system to help extend your hearing over a distance and in noisy environments.



ReSound Micro Mic and Multi Mic sits on the table or is clipped onto a person's clothing to stream sound directly to your hearing aids to improve hearing in noisy situations.



ReSound Phone Clip+ gives you more freedom, comfort and control while streaming phone calls, music or sound from your Bluetooth-enabled phone.



ReSound TV Streamer 2 streams stereo sound directly to your hearing aids, so you can listen at the volume that is comfortable for you.




ReSound Remote Control 2 allows you to adjust the volume or mute your hearing aids, switch programs, and see all your settings at a glance.



ReSound Remote Control intuitive control over volume and program settings with large, easy to use buttons.

Tinnitus is the #1 service- related injury among veterans



Find relief from tinnitus with ReSound

Tinnitus cannot be cured. But it can be managed. If you suffer from tinnitus, ReSound hearing aids provide you with new, smart tools for finding relief.

What is tinnitus?

From ringing and chirping to whooshing and clicking, tinnitus takes many forms. But most people refer to tinnitus as 'ringing in the brain'. For some, this ringing comes and goes without causing much bother. But for others, it is a constant noise that has a negative impact on daily life.

Manage your tinnitus

Using ReSound hearing aids helps change the way you respond to tinnitus by diverting your attention away from it with soothing background sounds. And, with the ReSound Smart 3D app, you can adjust these sounds, or choose from a range of sounds from nature to find comfort from your tinnitus:

Tinnitus Tuner

Find a personalized setting that provides you with relief from your tinnitus. Using the tuner to control the pitch and modulation of the sound in order to fine-tune the sound therapy that works best for you.

Nature Sounds

Try out the six pre-installed nature-inspired sounds, including 'Shoreline' and 'Beach Surf', to find out which one works best for you.



The Relief™ app

Direct audio streaming from iPhone, iPad and iPod touch gives you direct access to all the sound therapy apps available on the App Store, including ReSound Relief app, a multisensory app that uses a combination of customizable soundscapes and guided relaxation exercises to help distract you from your tinnitus.



ReSound GN

At ReSound we believe life is best when it's a natural, effortless flow, so we take to heart anything that holds you back. This is why we have devoted the past 150 years to making life sound better, connecting you to a life without limitations. Our breakthrough innovations, inspired by our Organic Hearing philosophy, all have one core ingredient – you, because hearing and hearing care should be as individual as a fingerprint.

ReSound is part of the GN Group – pioneering great sound from world-leading ReSound hearing aids to Jabra office headsets and sports headphones. The GN Group was founded in 1869, employs over 5,000 people, and is listed on NASDAQ OMX Copenhagen.

 facebook.com/resoundhearing

 twitter.com/resoundus

 youtube.com/resoundus

Manufacturer according to FDA:

ReSound Government Services

8001 E Bloomington Freeway
Bloomington, MN 55420
USA
1-800-392-9932
resound.com/veterans

© 2022 GN Hearing A/S. All rights reserved. ReSound is a trademark of GN Hearing A/S. Apple, the Apple logo, iPhone, iPad, iPod touch, Apple Watch and FaceTime are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google Inc. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories. The Bluetooth word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc.



MK604950 Rev G 2022.03